



HAFFKINE BIO PHARMACEUTICAL CORPORATION LIMITED

Procurement Cell

(A Government of Maharashtra Undertaking)

Regd. Office : Acharya Donde Marg, Parel, Mumbai 400 012 (INDIA)

Phone No: 022- 24129320-23

Website: <http://www.vaccinehaffkine.com>

Managing Director :022-24150628

E-mail: procurementcell@vaccinehaffkine.com

Corrigendum.

No.: 3984 / E-3130/ Change in Technical Specification / HBPCL/ PC/ 2021-22.

Date: 28/03/2022.

To,

M/s. Schiller Healthcare India Pvt. Ltd.,

Advance House, 2nd Flr, Makwana Road,

Off. Andheri Kurla Road, Marol Naka,

Andheri (E), Mumbai 59.

E-Mail: bhavin@schillerindia.com

Subject : - Change in Technical Specification

Reference: - 1.Tender no. E-3130/ Multipara Monitor with Accessories

2. 6502/Haffkine/Procurement Cell/E-3130/Multipara Monitor with Accessories/2021-22.


As per above mentioned / Change in Technical Specification

Instead of	Read As
1.The touch screen provides easy and intuitive operation. Just click the data area of the patient or parameter to see more information change a setting. Easy operation frees you to focus on your patient.	1.Multipara Monitor Should have touch screen and regulator knob. The touch screen provides easy and intuitive operation. Just click the data area of the patient or parameter to see more information change a setting. Easy operation frees you to focus on your patient.
2. Time is synchronized across all review screens so you can easily review an event at different points of view. Comprehensive review contributes to decision making for patients	2. Time is synchronized across all review screens so you can easily review an event at different points of view. Comprehensive review contributes to decision making for patients
3.With this interbed function, you can immediately check the status on any patient	3. With this interbed function, you can immediately check the status on any patient. CMS Connectivity Function is desirable.
4. This leads to better care of all patients.	4. This leads to better care of all patients.
5.With from any monitor you can check the information and alarm status on any other in the network	5.With from any monitor you can check the information and alarm status on any other in the network. CMS Connectivity function is desirable.
6.An onscreen guide with helpful illustrations makes it easy to operate the monitor and perform basic tasks even for the first time users	6.An onscreen guide with helpful illustrations makes it easy to operate the monitor and perform basic tasks even for the first time users. It is desirable that equipment demonstration CD and training to our staff on technical operational aspects of the equipment be provided at site, free of cost after installation and commissioning of the equipment at site.

7. Smooth operation enables you to focus on better patient care	7. Information brochure and guide should be provided.
8. Most effective parameter to prevent respiratory incidents	8. Deleted.
9. Early detection of trouble in breathing or ventilator	9. Most effective parameter to prevent respiratory incidents
10. Fast response with no lag time no warm up time compact and lightweight no heater or motor no sampling tube durable and easy to use for both non intubated and intubated patients	10. Monitors must have European CE (Notified body) or USFDA certification.
11. 10.4 inch color TFT LCD touchscreen	11. Early detection of trouble in breathing or ventilator.
12. Resolution: - 800×600 dots.	12. Fast Response with no lag time no warm up time compact and lightweight no heater or motor no sampling tube durable and easy to use for both non intubated and intubated patients. Should be durable, easy to use and fast response.
13. Number of waveforms: - ECG respiration, NIBP SPO ₂ temperature	13. Screen size should be atleast 12 inch color TFT LCD touchscreen.
14. Trend 120 hours.	14. Resolution: - 800×600 dots.
16. Numeric value list 120 hours	15. Number of waveforms: - ECG respiration rate, NIBP rate, SPO ₂ temperature, IBP
17. Arrhythmia recall 120 hours	16. Trend 120 hours.
18. Alarm history 120 hours.	17. 5 leads ECG, Pulse Oximeter, NIBP and Dual IBP
19. Battery operation time 6 hours.	18. Numeric value list 120 hours
20. Recorder 3 traces.	19. Should have Arrhythmia detection Function.
21. Network interface standard	20. Alarm history 120 hours.
22. Interbed 8 beds.	21. Battery operation time should be atleast 3-4 Hours.
23. Dimensions 330W×274H×156D mm	22. Recorder 3 traces
25. Weight 4 Kg.	23. Network interface standard
	24. Interbed 8 beds or 16 beds.
	25. Deleted
	26. Weight 5 kg or less

All other Terms & Conditions of the order will remain unchanged.

मा. व्यवस्थापकीय संचालक यांच्या मान्यतेने व करिता


Mrs. Sushama Patil
General Manager
(Additional Charge)

Haffkine Bio Pharmaceutical Corporation Ltd.
(Procurement Cell), Mumbai.